








Guangzhou Beulines Import & Export Co.,Ltd.

Add:No530,MinjieJiapin Commercial Building,Shatou Street,Panyu District,Guangzhou

Email:sales@beulines.com

Beufiller Cross-Linked Hyaluronic Acid Fillers

Types	Derm	Derm with lidociane	Deep	Deep with lidociane	Derm plus
Picture					
Volume	1ml /2ml	1ml	1ml/2ml	1ml /2ml	10ml / 20ml
Composition	24mg/ml hyaluronic acid				
Approximate number of gel particles 1ml	100,000		10,000		5,000
Particle size	0.15-0.28mm		0.28-0.5mm		0.5-1.25mm
Needle size	30G*2		27G*2		Not supplied
Where to inject	Middle part of dermis		Deep layer of dermis		Subcutaneous
Recommended indications	Glabellar lines,Tear trough,Nasolabial folds,Nasolabial groove,Oral commissures		Rhinoplasty,Nose,Lip enhancement,Chin augmentation,Cheeks padding,Cheek augmentation		Breast and Buttocks enlargement
Duration	6-12months		9-12 months		12-18 months
Storage	Room temperature($\leq 25^{\circ}\text{C}$), avoid to light.				
Origin	Non animal				
Color	Transport gel				
Package	One syringe and two needles of BD in a PVC blister, One blister and one instruction in a box.				

The price is different according to your order quantity .
PLS cousult us soon :)