



Guangzhou Beulines Import & Export Co.,Ltd.

Add:No530,MinjieJiapin Commercial Building,Shatou Street,Panyu District,Guangzhou
Email:sales@beulines.com

Beulines Cross-Linked Hyaluronic Acid Fillers

Types	Fine	Derm	Deep	Subskin
Volume	1ml / 2ml / 10ml / 20ml			
Composition	24mg/ml hyaluronic acid			
Approximate number of gel particles 1ml	100,000	100,000	10,000	5,000
Particle size	0.1-0.15mm	0.15-0.28mm	0.28-0.5mm	0.5-1.25mm
Needle size	30G*2	30G*2	27G*2	23G*2
Where to inject	Upper part of dermis	Middle part of dermis	Deep layer of dermis	Subcutaneous
Recommended indications	Fine lines Frown lines Crow's feet Lip lines	Deep wrinkles Nasolabial folds	Lip/Nose/Cheek/ Chin enlargement	Breast and Buttocks enlargement
Duration	6-12months		9-12 months	12-18 months
Storage	Room temperature($\leq 25^{\circ}\text{C}$), avoid to light.			
Origin	Non animal			
Color	Transparent gel			

**The price is different according to your order quantity .
PLS consult us soon :)**